

**Title:****Federated Skin Cancer Detection Using CNNs and FedOPT Optimization: A Privacy-Preserving Dermatology AI Framework with Client-Level Metric Tracking and 3D Learning Dynamics**

---

## Executive Summary

This paper presents a federated learning (FL) framework for skin cancer classification using the HAM10000 dermatology dataset. The system employs a lightweight convolutional neural network (CNN), integrates Flower's FedOPT optimization strategy, and incorporates client-level upsampling to mitigate data imbalance. The study also introduces comprehensive metric tracking and 3D visualization techniques to analyze per-client and per-round learning behavior. Results demonstrate stable convergence, improved classification accuracy under non-IID data conditions, and meaningful insights into federated optimization dynamics.

---

## 1. Introduction

Skin cancer is a major global health concern, and early detection significantly improves patient outcomes. Machine learning systems trained on dermoscopic images have shown promise in assisting clinical decision-making. However, the sensitive nature of medical data and regulatory constraints make centralized data aggregation impractical.

Federated Learning (FL) offers a decentralized alternative, allowing institutions to collaboratively train models without sharing raw patient data. This work develops an end-to-end FL pipeline for dermatology image classification, combining CNN modeling, adaptive federated optimization via FedOPT, client-level data balancing, and advanced evaluation through ROC curves and 3D surface visualizations.

---

## 2. Dataset: HAM10000

The HAM10000 dataset contains **10,015 dermoscopic images** across seven diagnostic categories:

- mel (melanoma)
- nv (melanocytic nevi)
- bkl (benign keratosis)
- df (dermatofibroma)
- akiec (actinic keratoses/Bowen's disease)
- bcc (basal cell carcinoma)
- vasc (vascular lesions)

### Challenges

- Severe class imbalance
- High intra-class variation
- Variability in image resolution and illumination

### Preprocessing

- Images resized to 128×128
- RGB normalization
- Label encoding
- Division into training and test sets
- Subsequent partitioning into five federated clients

Representative sample images illustrate the diverse visual characteristics across lesion categories.

---

## 3. Federated System Architecture

The federated learning system consists of five simulated clients and a central server orchestrated using the Flower framework.

### Key Components

#### 1. Client Nodes

- Receive unique, non-IID partitions of the training dataset
- Maintain separate training and validation splits
- Apply oversampling to balance class distribution

#### 2. Federated Server

- Performs FedOPT aggregation
- Tracks per-round metrics
- Updates and broadcasts global model parameters

#### 3. Training Protocol

- Five communication rounds
- Full participation from all clients each round
- Three local epochs per client

This architecture reflects real-world FL settings involving multiple healthcare institutions.

---

## 4. Model Architecture

A lightweight convolutional neural network (CNN) is employed to ensure computational efficiency in federated environments.

Conv2d(3, 32) → ReLU → MaxPool  
Conv2d(32, 64) → ReLU → MaxPool  
Conv2d(64, 128) → ReLU → MaxPool  
Flatten → Linear(32768, 128) → ReLU  
Linear(128, 7)

CrossEntropyLoss is used for multi-class classification, and the Adam optimizer is applied during local training.

---

## 5. Federated Optimization: FedOPT

FedOPT generalizes adaptive optimization techniques such as Adam to the federated setting.

### Parameters Used

- Server learning rate  $\eta = 0.01$
- Local learning rate  $\eta_l = 0.001$
- $\beta_1 = 0.90$
- $\beta_2 = 0.999$
- $\tau = 1e-9$

A customized FedOPT implementation is used to collect:

- Per-round training loss (per client)
- Per-round validation loss (per client)
- Per-round validation accuracy (per client)

These metrics enable fine-grained analysis of convergence behavior across heterogeneous clients.

---

## 6. Handling Non-IID Data Through Upsampling

The HAM10000 dataset exhibits substantial class imbalance. To reduce client drift and stabilize convergence, each client applies an oversampling procedure that:

- Identifies minority class samples
- Replicates them until all classes reach the majority-class frequency
- Randomizes dataset order after resampling

This strategy ensures more uniform class distribution across clients, improving both learning stability and evaluation consistency.

---

## 7. Training and Evaluation Methodology

### Local Training

Each client conducts:

- 3 training epochs
- Batch size of 32
- CNN forward and backward propagation
- Loss computation and optimization

### Local Validation

Each client computes:

- Validation accuracy
- Validation loss

## Global Evaluation

After federated training concludes, the final global model is evaluated on a held-out test set using:

- Confusion matrix
  - Class-wise ROC curves
  - AUC scores
- 

# 8. Results


## 8.1 Accuracy Surface

A 3D surface visualization of validation accuracy across rounds and clients shows:

- Consistent improvement over all five rounds
- Final-round accuracies ranging between  $\sim 0.80$  and  $\sim 0.84$
- Reduced variability among clients in later rounds

### Interpretation:

FedOPT demonstrates stable convergence even with heterogeneous client datasets.



(Figure: Accuracy Surface over Rounds & Clients)


---

## 8.2 Loss Surface

Training and validation losses decrease steadily across rounds, with final validation losses in the 0.35–0.45 range.

### Interpretation:

The downward trend indicates effective learning dynamics with no signs of divergence or overfitting.



(Figure: Loss Surface over Rounds & Clients)

## 8.3 Multi-Class ROC and AUC


One-vs-rest ROC analysis yields the following AUC values:

Classes	AUC
df	0.78
mel	0.57
bkl	0.55
bcc	0.51
nv	0.44
vasc	0.39
akiec	0.32

Interpretation:



- Strong separability for *df*
- Moderate performance for *mel* and *bkl*
- Lower performance for *akiec* and *nv* due to high intra-class variability and limited samples



(Figure: Multi-class ROC Curve)

## 8.4 Sample Images

A grid of representative images demonstrates the inherent complexity of the classification task across lesion types.



*(Figure: Sample Images from Dataset)*

---

## 9. Discussion

The results demonstrate that:

- The FedOPT-based FL system converges reliably on non-IID medical image data
- Upsampling is effective in reducing performance disparities among clients
- 3D visualization is highly useful for diagnosing cross-client learning dynamics
- Despite architectural simplicity, the CNN achieves competitive performance

## Limitations

- Minority classes remain challenging due to limited sample sizes
  - A simple CNN architecture constrains maximum achievable accuracy
  - Further improvements require more sophisticated model architectures or augmentation techniques
- 

## 10. Future Work

Potential extensions include:

- Differential privacy integration for regulatory compliance
  - Secure aggregation for enhanced confidentiality
  - Personalized federated methods (e.g., FedPer, FedRep)
  - Vision Transformers adapted for federated settings
  - Automated hyperparameter tuning
  - Deployment across real multi-institution networks
- 

## 11. Conclusion

This study presents a practical and interpretable federated learning pipeline for skin cancer detection using non-IID dermatology image data.

The combination of:

- Lightweight CNN modeling
- FedOPT optimization

- Client-level upsampling
- Per-client metric tracking
- 3D convergence visualizations

collectively demonstrates the feasibility and value of applying FL to medical imaging tasks.

The framework serves as a foundation for further research and development in privacy-preserving clinical AI systems.

---

## 12. References

- Flower: A Friendly Federated Learning Framework
  - Tschandl et al. HAM10000 Dataset
  - Reddi et al., Adaptive Federated Optimization Algorithms (FedOPT)
  - PyTorch Documentation
  - scikit-learn Metrics Documentation
-